NISSEI BLOOD PRESSURE MONITOR WSK-1011

INSTRUCTIONS

This manual is intended to assist you in the safe and efficient operation of BLOOD PRESSURE MONITOR WSK-1011. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

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
PL

CONTENTS

| | |
|------------------------------------------|---|
| PART NAMES AND PRODUCT COMPONENTS . . . | 2 |
| GENERAL INFORMATION | 2 |
| PRECAUTIONS AND CONTRAINDICATIONS . . . | 5 |
| TIPS ON TAKING YOUR BLOOD PRESSURE . . . | 6 |
| PREPARATION BEFORE USE | 7 |
| MAKING A MEASUREMENT | 9 |

| | |
|--------------------------------------------------------|----|
| VIEWING THE SAVED READINGS (MEMORY FUNCTION) | 12 |
| TROUBLESHOOTING | 13 |
| TECHNICAL DESCRIPTION | 14 |
| WARRANTY | 96 |

PART NAMES AND PRODUCT COMPONENTS

CARRYING BAG

Use the bag to store the monitor to avoid accidental activation when carrying.

AAA (LR03) BATTERIES

For test use

INSTRUCTIONS

GENERAL INFORMATION

INDICATIONS FOR USE

This product is intended for noninvasive measurement of systolic and diastolic

blood pressure, determination of pulse rate and calculation of pulse pressure in adults in a home healthcare environment. The product is not designed for

neonatal use. Please consult with your doctor or physician to use this product to take blood pressure of child or person in pregnancy or under pre-eclamptic condition.


METHOD OF MEASUREMENT

This product employs the oscillometric method for measurement of blood pressure and pulse rate. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the wrist in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value.

FUZZY INFLATION

This product employs FUZZY inflation to ease stress on wrist during cuff inflation.

By detecting pulse signals while inflating, this product is able to determine the sufficient pressure for each measurement.

WHO CLASSIFICATION SYMBOL

Measured blood pressure values are classified against WHO guideline (1999). The higher class applies when systolic and diastolic pressures fall into different classes.

For specific information on your blood pressure, contact your physician. Never make any judgment on your own based on the classification indication only.

| DISPLAY | WHO CLASSIFICATION | SYS | DIA |
|---------|-------------------------|---------|---------|
| | Hypertension (severe) | ≥180 | ≥110 |
| | Hypertension (moderate) | 160-179 | 100-109 |
| | Hypertension (mild) | 140-159 | 90-99 |
| | High Normal | 130-139 | 85-89 |
| | Normal | 120-129 | 80-84 |
| | Optimal | <120 | <80 |

PULSE PRESSURE DISPLAY

This product calculates and displays pulse pressure.

While systolic blood pressure keeps increasing with age, diastolic blood pressure tends to start decreasing at around age 50. Pulse pressure is value subtracting diastolic value from diastolic value and therefore pulse pressure tends to increase with age. Pulse pressure has been considered to be related to stiffness of artery and has been studied as one of circulatory risk factors. It is sometimes said that 45 mmHg of pulse pressure is normal value. Pulse pressure only does not represent arterial sclerosis. However, observation of pulse pressure in long term will be important.


BODY MOTION SYMBOL

Blood pressure value taken while moving cannot be said to be the correct value because body movement can affect blood pressure.

This product analyzes pulse wave and displays [] when body motion is detected. [] indicates the results might be affected by body movement.

IRREGULAR PULSE RHYTHM SYMBOL

Pulse rhythm can be disturbed from talking, moving or arrhythmias. This product displays [] when the difference of shortest and longest time intervals between pulse waves is over 25%, indicating irregular pulse rhythm.

MEMORY FUNCTION

The measured values are automatically saved for later review in either one of two memory banks, [1] or [2]. These two memory banks can be used to save readings of two persons separately or to save morning readings and evening readings separately.

Each bank can save up to 60 readings. When the number of stored readings reaches 60, the oldest reading will be deleted to record new reading.

Saved readings are recalled with measurement date and time when clock function is activated. *The clock needs to be activated before measurement for date and time to be saved with the measured valued.

[Err] results are not saved.

CARE AND MAINTENANCE

Keep the product clean. Inspect its cleanliness after use. To clean, use only a soft dry cloth. Do not use gasoline, paint thinner, or other strong solvents. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, use a synthetic detergent and gently rub the surface.

When storing the product, do not place heavy objects on it. Do not fold or bend the cuff forcibly. When the product has been stored at a temperature below the freezing point, keep it for at least 1 hour in a warm place before using. Remove the batteries if the product is to be stored for an extended period of time. Keep the batteries out of reach of children.

We suggest that you have your blood pressure monitor checked every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

PRODUCT SPECIFICATIONS

| | |
|---------------------------|-------------------------------------------------------------------------------------|
| Model | : WSK-1011 |
| Operating Principle | : Oscillometric method |
| Indicator | : 15 digits liquid crystal display |
| Pressure Indicating Range | : 3 to 300 mmHg (cuff pressure) |
| Measuring Range | : 50 to 250 mmHg (systolic), 40 to 180 mmHg (diastolic), 40 to 160 bpm (pulse rate) |
| Accuracy* | : ± 3 mmHg (cuff pressure), $\pm 5\%$ of reading (pulse rate) |
| Inflation | : Automatic inflation (FUZZY logic) |
| Deflation | : Automatic (electric control valve) |

| | |
|---------------------------------------------|----------------------------------------------------------------|
| Exhaust | : Automatic exhaust valve |
| Power Supply | : Two 1.5 volt AAA alkaline (LR03) type batteries |
| Power Consumption | : 2W (max.) |
| Electric Rating | : DC3V/2W |
| Memory | : 2 banks, each stores 60 readings |
| Operating Condition | : +10°C to +40°C, 30% to 85% RH (noncondensing) |
| Transportation/Storage Condition | : -20°C to +60°C, 10% to 95% RH (noncondensing) |
| Coverage wrist circumference | : 12.5 to 22.5 cm |
| Weight | : Approx. 110 g (without batteries) |
| Size | : Approx. 70 x 27 x 70 mm (W x D x H) , not including the cuff |
| Protection against harmful ingress of water | : IPX0 (IEC 60601-1) |
| Protection against electric shock | : Internally powered equipment, Type BF applied part (cuff) |
| Mode of operation | : Continuous operation |

| | |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Classification | : Internally powered equipment |
| Key to symbols | :  : Type BF applied part |
| | :  : Refer to instruction manual/booklet. |
| | :  : The used electrical and electronic products are not household waste. Follow your national/local recycling rules to dispose of them properly. In the EU countries, please refer to waste management symbol(s) marked on the package or the instrument. |

This device complies with EN1060-1:1995+A2:2009 Non-invasive sphygmomanometers Part 1: General requirements, EN1060-3:1997+A2:2009 Non-invasive sphygmomanometers Part 3: Supplementary requirements for electro-mechanical blood pressure measuring system and EMC (IEC60601-1-2:2007).

*Accuracy is guaranteed with the measured values that are within the measuring range.

Specifications are subject to change without notice due to improvements in performance.

PRECAUTIONS AND CONTRAINDICATIONS

Do not use this product without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.

For specific information on your blood pressure, contact your

physician. Never make any judgment on your own regarding measurement results.

Do not use this product along with implantable and wearable medical electrical equipment such as pacemaker, defibrillator, or electrocardiographic monitor. This product is not also intended to be

used with HF surgical equipment.

Do not use this product in an explosive environment such as near flammable anesthetics or inside oxygen chamber.

The system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.

Do not use any accessories other than those specified by the manufacturer. Otherwise, correct measurement readings cannot be obtained.

Use of this product in areas near mobile phones, microwave ovens or other devices with strong electromagnetic field may cause malfunctions.

Do not apply the cuff over wounded wrist, wrist of which side is under an intravascular access or therapy or an arterio-venous shunt, or wrist on the side of a mastectomy. Otherwise injury may be resulted.

Make sure that inflation of the cuff is not causing prolonged

impairment of blood circulation. Also, be cautious about temporary loss of the functions of any other medical equipment if any monitoring equipment is used on the same limb with the blood pressure measuring cuff.

To avoid harmful injury due to interfered blood flow from cuff inflation, do not make measurements repeatedly.

Because the product includes precision parts, avoid extreme temperature variations, humidity, shock, dust, and direct sunlight. Do not drop or strike the product. Make sure not to expose it to moisture. This product is not water resistant.

Do not press the display and operation keys or place the monitor with display face down.

Do not take out batteries when the monitor is turned on. Make sure to switch off the monitor before removing batteries.

Do not disassemble or modify the product.

Do not inflate the cuff when it is not wrapped around your wrist.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins in this flow.

Blood pressure varies all the time, influenced by mental and physical factors, and is never constant.

In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and increases at a relatively fast rate after arising from bed.

Causes for Changes in Blood Pressure

- Body movement
- Mental Tension
- Emotions
- Drinking Alcohol
- Recent Urination or Bowel Movement
- Changes in the surroundings such as movement or noise, etc.
- Conversation
- Nervousness
- Eating
- Smoking
- Room temperature

Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.

This is because you are tense at the hospital and relaxed at home. It is important to know your stable normal blood pressure at home.

Blood pressure at the wrist may be different from that at the upper arm.

As the thickness of the blood vessels differs depending on measuring position, blood pressure value measured at the wrist may be slightly different from that measured at the upper arm. In most cases, the difference for a healthy person will be within about ± 10 mmHg for the systolic and the diastolic blood pressure. Caution is required for the following persons: persons with hypertension, diabetes, liver trouble, hardened arteries, peripheral hindrance of circulation, etc., as there may be a difference between the values measured at the wrist and those measured at the upper arm.

Let a qualified physician interpret your blood pressure readings.

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

Before taking blood pressure, rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.

Measure blood pressure using the correct posture and do not move or speak during measurement.

Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement.

Take your blood pressure at the same time every day.

The ambient temperature should be approximately 20°C when you take your blood pressure.

Reading may vary slightly, depending on posture, whether sitting or lying.


PREPARATION BEFORE USE

Install batteries and activate the clock before use. Measurement can be made with clock being inactivated. However, activate the clock to

save measurement date and time with readings for better review. The clock is inactivated as the batteries are removed.

1 INSTALLING BATTERIES

Open the battery compartment by unlocking the cover by pressing its top and bottom.

Install two identical AAA alkaline (LR03) batteries.

Make sure that the polarities correspond to the (+) and (-) marks inside the battery compartment.

Batteries can be easily installed or removed as their (-) ends are pushed against the spring.


You may use nickel hydride rechargeable batteries.

NOTE: For details regarding use of rechargeable batteries, follow the directions accompanying them.

You will hear a beep as the batteries are installed.

Close the battery compartment cover.

Do not force the battery cover into position.

**AAA (LR03)
BATTERIES**


2 ADJUST DATE AND TIME AND ACTIVATE THE CLOCK

Touch "SET" KEY and keep your finger on the key until "20 11" starts flashing.

Clock is set in the order of year, month, day, hour and minute.

The flashing number increases with MEMORY KEY [1] and decreases with KEY [2]. The number will be fast-forwarded if you keep your finger on the key. Touching "SET" key will fix the number and the next article will flash. Touching "START/STOP" key will terminate the setting.

As the date and time are set, clock is activated. Display of the clock while the device is turned off indicates that the clock is activated.

Battery symbol [] will flash when battery power is low. Measurement can not be made while the symbol is displayed. Replace the batteries. Both batteries need to be replaced with new ones; do not mix new and old batteries or different batteries.

The symbol might appear only during measurement although you can review readings in memory. This is because more battery power is required for conducting a measurement, e.g. inflating the cuff, than displaying the

readings in memory. Please have the batteries replaced.
The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries.

MAKING A MEASUREMENT

The cuff is applicable to wrist with circumference between 12.5 and 22.5 cm. Take your wrist size before measurement.


1 SIT ON A CHAIR AND APPLY CUFF.

Find a chair and a table so that you can be comfortably seated with your feet flat on the floor and your back and arm supported.

Remove watch or bracelet.

Put the cuff on your left wrist directly on the skin with the display screen on the side of your palm.

The monitor should be placed at the center of your wrist with the top of the monitor positioned 5 to 10 mm from the bottom of your palm.


Fasten the cuff and press the surface of the cuff to make sure the hook and loop attaches securely.

The cuff should fit snugly and there should be no space between the cuff and your wrist.


Take care that your clothing is not caught beneath the cuff.

You may use your right wrist for measurement. However, make sure to use same wrist for measurement; the blood pressure on left wrist and right wrist may not be the same.

2 ADJUST THE MONITOR TO THE HEIGHT OF YOUR HEART.

Place a folded towel or similar under your arm to support it so that the cuff is kept at the height of your heart during measurement.

Hold your elbow if you do not have an appropriate table to put your elbow on.

Measurement can be made lying on your back. Make sure that the cuff is at the height of your heart.


If the cuff is lower (higher) than the heart, the measured reading tends to become larger (smaller).

Take deep breaths and relax.

Do not move or speak while taking your blood pressure.

Do not cross your legs during measurement.

Do not hold your breath.

3 TURN ON THE MONITOR.

Touch “START/STOP” KEY.

Air is exhausted from the cuff.


Automatic inflation starts.

Inflation stops at optimal pressure and pressure starts to decrease.


PULSE SYMBOL is displayed as pulse is detected.

Touch “START/STOP” KEY again to cancel measurement.

The monitor will exhaust air from the cuff and turn off.

Air is released from the cuff and measured values are displayed as measurement completes.

Refer to page 4 for [].

Refer to page 4 for [].

Refer to page 13 for “Err” display.

AUTOMATIC REPEATED INFLATION

If you move or strain your arm or hand, inflation may be repeated multiple times to obtain a measurement. Repeated inflation does not indicate malfunction of the monitor.

4 SELECT MEMORY BANK TO SAVE THE MEASUREMENT RESULT.

MEMORY SYMBOL is displayed below MEMORY KEY.

Touching **MEMORY KEY** will switch the bank.

MEMORY MARK

Readings are automatically saved in the bank, either MEMORY [1] or [2], selected when the monitor is turned off after measurement. Make sure to confirm the bank selected after a measurement. Refer to page 4 and page 12 for memory function.

5 TURN OFF THE MONITOR.

Touch “**START/STOP**” KEY.

The monitor will be turned off. Even if you do not turn off the monitor, it turns off automatically after 3 minutes.

Do not execute repeated measurements. This will cause engorgement of the wrist and correct values will not be obtained. Let your wrist rest for at least 5 minutes.

VIEWING THE SAVED READINGS (MEMORY FUNCTION)


Turn off the monitor once with "START/STOP" KEY when the results are displayed after a measurement.
Refer to page 4 for memory function.

1 TOUCH MEMORY KEY.

Touch MEMORY KEY [1] to review readings saved in the memory bank [1] and MEMORY KEY [2] to review those in the bank [2].

The average of the saved readings are displayed.


Note: Average will not be displayed unless there are two or more readings saved.

2 TOUCH THE SAME MEMORY KEY AGAIN.

The latest reading saved in the selected bank is displayed.

The display at the top shows memory number, measurement date and time alternately. The latest reading is displayed as memory "No. 1".

The results are recalled from the latest to the oldest with operation of MEMORY KEY.

The bigger memory number indicates the older number.

SWITCHING MEMORY BANKS

Touching MEMORY KEY [2] at display of reading in the bank [1] will switch the display to reading in bank [2] and touching MEMORY KEY [1] again will return to the display of result in the bank [1].

3 TURN OFF THE MONITOR.

Touch "START/STOP" KEY.

Even if you do not turn off the monitor, it displays the result for approximately 30 seconds and will turn off automatically.

DELETING SAVED READINGS

1 DISPLAY THE READING TO DELETE.

Display the average reading to delete entire readings in the memory.

2 TOUCH MEMORY KEY AND KEEP YOUR FINGER TOUCHED TO THE KEY.

Touch MEMORY KEY [1] to delete reading in the memory bank [1] or MEMORY KEY [2] to delete reading in the bank [2].

The display will start flashing. Keep your finger on the key until the reading is replaced with " - - -".

TROUBLESHOOTING

???

Blood pressure is extremely high or low.

The monitor is not adjusted to your heart. **Adjust the monitor to the height of your heart.**

Measurement is not taken with correct posture or the cuff is wrapped incorrectly. **Reconfirm the measurement procedures.**

The reading was effected from moving or talking. **Remain still and quiet during measurement.**

Measurement was taken right after exercise or other activities influencing the reading. **Measure again after resting for more than 5 minutes. See TIPS ON TAKING YOUR BLOOD PRESSURE.**

???

Measured values vary all the time.

Mental and physical condition or measurement condition is influencing blood pressure. **Make measurement under the same condition.**

???

Measured value is different from that taken at a hospital.

Mental condition such as nervousness at a hospital or relaxed feeling at home is influencing blood pressure. **Take records of blood pressure taken at home and consult with your doctor.**

Blood pressure is changing 24 hours a day in relation with mental and physical conditions. Eating, drinking, smoking, exercising or bathing will alter blood pressure. Blood pressure will be also influenced by mental condition such as feeling tensed or relieved and health condition It is important to know the tendency of blood pressure by making measurement regularly. Make your own rules including conditions and time to make measurement and take blood pressure everyday.

???

Inflation is repeated.

The cuff is inflated again if sudden change of pressure, which could be due to movement or straining of arm or wrist, was detected. Inflation is repeated until blood pressure is taken. **Repeated inflation is not a malfunction of the monitor. Do not move or speak during measurement.**

Err

300

is displayed.

Over pressure: blood pressure could not be taken due to moving or speaking although the cuff was inflated to the maximum pressure. **Do not move or speak during measurement.**

Err

-1

or

Err

-3

is displayed.

Blood pressure could not be taken due to moving or speaking. **Do not move or speak during measurement.**

is displayed.

The cuff is not properly applied. **Properly apply the cuff. If Err-2 appears after re-trying to apply the cuff properly, assume the cuff is damaged. Contact your dealer.**

is displayed.

Batteries are weak. **Replace all batteries with new ones.**

Nothing is displayed.

Batteries are depleted. **Replace all batteries with new ones.**

Batteries are inserted incorrectly. **Reinsert the batteries correctly.**

Battery terminals are not clean. **Clean the terminals with dry cloth.**

Touching the keys with tip of finger or nail is touching the key. **Touch the key with the plane of your finger.**

Clock is not displayed.

The clock is inactivated. NOTE: Clock is inactivated when batteries are removed. **Adjust date and time and activate the clock.**

Measurement date and time are displayed with "--/--" and "--:--".

The clock is not activated. **Adjust date and time and activate the clock.**

Measurement was taken before the clock was activated. **Measurement date and time cannot be saved without clock being activated.**

is displayed.

START/STOP KEY was being touched at installation of batteries. **Turn off the monitor once with "START/STOP" KEY.**

If you cannot get correct measurement with the methods above, contact your dealer. Do not disassemble or tamper with the internal mechanism.

TECHNICAL DESCRIPTION

WSK-1011 complies with the EMC, electromagnetic compatibility, standard, IEC60601-1-2. Refer to the following tables for specific information regarding compliance to the standard. WSK-1011, as a medical electrical equipment, needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

- Portable and mobile RF communications equipments can affect the device.
- The use of accessories other than those specified in this manual may result in increased emissions or decreased immunity of the device.
- WSK-1011 should not be used adjacent to or stacked with other equipment.

Guidance and manufacturer's declaration - electromagnetic emissions -

| WSK-1011 is intended for use in the electromagnetic environment specified below. The customer or the user of WSK-1011 should assure that it is used in such an environment. | | |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Emissions test | Compliance | Electromagnetic environment - guidance |
| RF emissions CISPR 11 | Group 1 | WSK-1011 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |

| | | |
|------------------------------------------------------|---------|---------------------------|
| RF emissions CISPR 11 | Class B | Note: Operate by DC Power |
| Harmonic emissions IEC 61000-3-2 | N/A | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | N/A | |

Guidance and manufacturer's declaration - electromagnetic immunity -

WSK-1011 is intended for use in the electromagnetic environment specified below. The customer or the user of WSK-1011 should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|-----------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines | N/A | Note: Operate by DC Power. |
| Surge IEC 61000-4-5 | ±1 kV differential mode ±2 kV common mode | N/A | Note: Operate by DC Power. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_n (>95% dip in U_n) for 0,5 cycle 40% U_n (60% dip in U_n) for 5 cycles 70% U_n (30% dip in U_n) for 25 cycles <5% U_n (>95% dip in U_n) for 5 sec | N/A | Note: Operate by DC Power. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

NOTE U_n is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity -

WSK-1011 is intended for use in the electromagnetic environment specified below. The customer or the user of WSK-1011 should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|---------------|----------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| | | | Portable and mobile RF communications equipment should be used no closer to any part of WSK-1011, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. |

| | | | |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 3 Vrms, 150 kHz to 80 MHz 3 V/m, 80 MHz to 2.5 GHz | 3 V 3 V/m | Recommended separation distance $d=[3.5/V1] \sqrt{P}$ $d=[3.5/E1] \sqrt{P}$, 80MHz-800MHz $d=[7/E1] \sqrt{P}$, 800MHz-2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:  |
| NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | |
| a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radios broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which WSK-1011 is used exceeds the applicable RF compliance level above, WSK-1011 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating WSK-1011. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than $[V1] \sqrt{P}$. | | | |

Recommended separation distances between portable and mobile RF communications equipment and WSK-1011 -

WSK-1011 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of WSK-1011 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and WSK-1011 as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter, W | Separation distance according to frequency of transmitter, m | | |
|----------------------------------------------|--------------------------------------------------------------|---------------------------------------------|--------------------------------------------|
| | 150 kHz to 80 MHz, $d=[3.5/V1] \sqrt{P}$ | 80 MHz to 800 MHz, $d=[3.5/E1] \sqrt{P}$ | 800 MHz to 2.5 GHz, $d=[7/V1] \sqrt{P}$ |
| 0,01 | 0.12 | 0.12 | 0.23 |
| 0,1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

This equipment is guaranteed for the period of 2 years after the date of purchase against manufacturing defects when returned along with the proof of date of purchase to the dealer from whom the purchase was made. During this period, the unit will be repaired or replaced free of charge if the failure is attributable to faulty design or manufacture. This warranty does not cover damage or malfunctions caused by improper handling or use contrary to the instructions in this manual. Please contact your dealer for additional warranty provisions which may remain effective after the manufacturer's warranty period has expired.

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LTD.**
JPI Inc EU Office